

# Total Mercury in Soils

LECO Corporation; Saint Joseph, Michigan USA

## Instrument: AMA254



### Accessories

614-822-102 Small Nickel Boats

*Note: Boats should be pre-baked at 400°C or analyzed (without a sample) before loading a sample.*

### Sample Weight

25 to 80 mg (0.025 to 0.080 g)

### Calibration Standard

LECO 502-813 Fly Ash, LECO 502-499 (BCR 143r), LECO 502-649 Dry Sludge (NIST 2781), or other suitable reference material

### Analysis Time

~8 minutes

### Method Profile

Drying Time:	60 seconds
Decomposition Time:	200 seconds
Cuvette Clear Time:	45 seconds
Dosing Delay Time:	0 seconds
Cell Selection:	Auto select
Metric for Calculations:	Peak Area

*NOTE: Method for Quicksilver Windows® Software Version 2.0.*

### Procedure

- Determine the blank as follows.
  - Enter "Blank" from the drop-down menu under the "Name" column.
  - Click "Analyze"; the door will open and the nickel loop will be presented.
  - Carefully place a 614-822-102 Small Nickel Boat into the nickel loop using clean tweezers.
  - Click "OK" in the "Load Sample" window; the door will close and the analysis sequence will start automatically.
  - Repeat steps 1a through 1c two more times. The system and boats will be purged of any interfering elements.
- Calibrate the instrument as defined in the instructional manual.
  - Weigh various weights in accordance to the absolute amount of mercury required to calibrate an appropriate dynamic range. The calibration samples are weighted into the 614-822-102 Small Nickel Boat.
  - Enter each calibration sample with their appropriate ID code from the drop-down menu, and sample weight from an external balance measurement.
  - Click "Analyze"; the door will open and the nickel loop will be presented.
  - If there is a boat in the nickel loop, remove it and keep for later use.
  - Carefully place the calibration sample boat into the nickel loop using clean tweezers.
  - Click "OK" in the "Load Sample" window; the door will close and the analysis sequence will start automatically.
  - Repeat steps 2a through 2f as per the calibration procedures.

*Note: The first analyzed sample after a long delay should be discarded. This sample should be considered a conditioner for the system, and not used for the actual calibration.*
- Complete a calibration by following the calibration procedure as outlined in the manual.
  - Verify the calibration by analyzing one of the calibration samples again. It should be within the expected tolerances. If not, repeat steps 2a through 2i again.
- Analyze the samples as follows.
  - Weigh ~50 mg of the sample into a 614-822-102 Small Nickel Boat.
  - Enter a sample identification in the Name column and the sample weight in the Mass column.
  - Click "Analyze"; the door will open and the nickel loop will be presented.
  - If there is a boat in the nickel loop, remove it and keep for later use.
  - Carefully place the sample boat into the nickel loop using clean tweezers.
  - Click "OK" in the "Load Sample" window; the door will close and the analysis sequence will start automatically.

### Typical Results

Sample	Weight(mg)	Hg (ng)	Hg (ppm)
NIST 2709	47.3	67.9	1.435
Soil (as is)	50.5	72.5	1.435
	48.1	72.1	1.479
	48.3	69.1	1.432
	49.5	71.1	1.436
		<b>Avg (ppm)</b>	<b>1.44</b>
		<b>Std</b>	<b>0.02</b>
		<b>RSD</b>	<b>1.39%</b>

Sample	Weight(mg)	Hg (ng)	Hg (ppm)
Soil #1	76.8	17.1	0.223
	77.5	16.6	0.214
	75.9	16.5	0.217
	78.4	17.6	0.224
	74.3	16.5	0.222
		<b>Avg (ppm)</b>	<b>0.220</b>
		<b>Std</b>	<b>0.004</b>
		<b>RSD</b>	<b>1.96%</b>

Sample	Weight(mg)	Hg (ng)	Hg (ppm)
Soil #2	23.2	166.5	7.18
	25.8	174.0	6.74
	23.7	171.7	7.25
	22.4	159.9	7.14
	29.2	194.6	6.70
		<b>Avg (ppm)</b>	<b>7.00</b>
		<b>Std</b>	<b>0.261</b>
		<b>RSD</b>	<b>3.7%</b>



### LECO Corporation

3000 Lakeview Avenue • St. Joseph, MI 49085 • Phone: 800-292-6141 • Fax: 269-982-8977  
 info@leco.com • www.leco.com • ISO-9001:2000 • No. FM 24045 • LECO is a registered trademark of LECO Corporation.